



70° ALCOHOL RANGE

Dermatologically tested disinfectant spray for external use.

- ✓ **Body disinfectant (on healthy skin)**
- ✓ **Surface disinfectant (small equipment)**

DISINFECTANT PROPERTIES IN COMPLIANCE WITH STANDARDS AT 20°C:

| Standards | Strains | Conditions |
|--|---|--------------|
| Bactericidal in compliance with EN 1276 | <i>Staphylococcus aureus</i> , <i>Escherichia coli</i> , <i>Pseudomonas aeruginosa</i> , <i>Enterococcus hirae</i> | In 5 minutes |
| Bactericidal in compliance with EN 13727 | | In 1 minute |

APPLICATION

As antiseptic on healthy skin.
Can also be used on small equipment: tweezers, small scissors, etc.

INSTRUCTIONS FOR USE

Apply disinfectant to healthy skin. Avoid contact with eyes.
Leave for 5 minutes.

PACKAGING

A13832 - 70° ALCOHOL 2 L

Presentation: 2 L container
Pack size: Box of 4 units
Pallet size: 300 units – 75 boxes

A13831 - 70° ALCOHOL 1 L

Presentation: 1 L container
Pack size: Box of 12 units
Pallet size: 600 units – 50 boxes

A12831 - 70° ALCOHOL 50 mL

Presentation: 50 mL aerosol
Pack size: Box of 24 units
Pallet size: 4608 units – 192 boxes

REGULATIONS

FOR INDIVIDUAL AND PROFESSIONAL USE
Biocide category: TP1, human hygiene.
Dangerous. Follow precautions for use.
Use biocides with precaution. Read the label and other information about the product before use.
Our data sheets are prepared with current knowledge and subject to modifications. Photo non-binding.
29/10/2024 Index 3.2



Retrouvez toutes nos fiches techniques sur www.sico.fr

SIEGE SOCIAL / HEADQUARTER
577 Rue du Pommarin B.P. 16 – Voreppe
38341 Moirans CEDEX France
Tél : +33 (0)4 76 50 85 50 / sico.commercial@sico.fr

